

EVENT MARKER/ WATER PROOF

Ref. BCTOCO16

COMEN



CE
0482



EVENT MARKER/ WATER PROOF

COMEN

Ref. [BCTOCO16](#)

C11 cardiocograph probe to monitor mom's contractions

Description

The event marker for the Comen C11 cardiocograph allows the patient to manually indicate perceived fetal movements during monitoring. It helps healthcare professionals record key clinical events during pregnancy and labor surveillance, supporting better data interpretation.

Features

CE Marking	CE0482	Medical device	Class Ila
Bar code	3661040219271	HS Code	9018193010
Warranty (months)	6	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

02/04/2026 23:43:08 – For information purposes only, can be subject to change