

## RING CUTTER

Ref. IP05420



02/04/2026 23:59:47 – For information purposes only, can be subject to change

# RING CUTTER

Ref. IP05420



**Used to easily cut a ring safely**

## Description

- Cut materials traditionally used in jewelry
- Gold, silver, brass, copper, with a thickness up to 2.5 mm

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

100% IN-STORE reference: this reference of service instrumentation is guaranteed to be available immediately.

## Features

CE Marking	NONE	Medical device	No
Bar code	3661040021935	HS Code	90189084
Warranty (months)	24	Country	Pakistan

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	23	Bare Unit / Width (cm)	0.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	6
Packed unit / Length (cm)	9	Bare Unit / Length (cm)	19
Packed unit / Weight	0.14	Bare Unit / Weight (kg)	0.136

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 02/04/2026 23:59:47 – For information purposes only, can be subject to change