

COTTON BOX 11X14 CM

Ref. AI01114



03/04/2026 00:42:38 – For information purposes only, can be subject to change

COTTON BOX 11X14 CM

Ref. [AI01114](#)



Allows cotton to be stored to protect against dust

Description

- Diameter 110 cm
-
- Sold with lid for effective and hygienic protection
-
- Decontamination with soapy water (do not use detergent) -
- Sterilizable in autoclave (134°-18 minutes) or in dry air
-
- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
-
- 100% IN-STORE reference: this service instrumentation reference is guaranteed to be available immediately

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019208	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	12	Bare Unit / Width (cm)	11
Packed unit / Height (cm)	12	Bare Unit / Height (cm)	14
Packed unit / Length (cm)	17	Bare Unit / Length (cm)	14
Packed unit / Weight	0.362	Bare Unit / Weight (kg)	0.362

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

 03/04/2026 00:42:38 – For information purposes only, can be subject to change