

**S/S BED PAN WITHOUT COVER SLIPPER
TYPE 48 X 30CM**

Ref. AI14830



Description

Capacity 6 Liters

Maximum weight supported: 90kg

03/04/2026 00:05:13 – For information purposes only, can be subject to change

S/S BED PAN WITHOUT COVER SLIPPER TYPE 48 X 30CM



Ref. AI14830

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019352	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	55	Bare Unit / Width (cm)	30
Packed unit / Height (cm)	8	Packed unit / Length (cm)	30
Bare Unit / Length (cm)	48	Packed unit / Weight	1.08
Bare Unit / Weight (kg)	1.07		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

03/04/2026 00:05:13 – For information purposes only, can be subject to change