

S/S BOX 20 X 10 X 3 CM

Ref. AI20103



03/04/2026 01:24:07 – For information purposes only, can be subject to change

S/S BOX 20 X 10 X 3 CM

Ref. AI20103



Convenient for storing surgical instruments

Description

- Sold with lid to protect instruments
- Decontamination with detergent and disinfectant -
- Sterilizable by autoclave (134°-18 minutes)- In the " ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- Reference 100% IN-STORE: this service instrumentation reference is guaranteed available immediately
- Available in different sizes

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019499	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	10	Bare Unit / Width (cm)	10
Packed unit / Height (cm)	3	Bare Unit / Height (cm)	3
Packed unit / Length (cm)	20	Bare Unit / Length (cm)	20
Packed unit / Weight	0.401	Bare Unit / Weight (kg)	0.4

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 03/04/2026 01:24:07 – For information purposes only, can be subject to change