

**S/S BOX PERFORATED 20 X 10 X 5CM
PERFORATED TOP / BOTTOM**

Ref. AI20105P



03/04/2026 00:00:07 – For information purposes only, can be subject to change

**S/S BOX PERFORATED 20 X 10 X 5CM
PERFORATED TOP / BOTTOM**



Ref. AI20105P

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040067742	HS Code	90189084
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	12	Bare Unit / Width (cm)	10.5
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	5
Packed unit / Length (cm)	21	Bare Unit / Length (cm)	20.7
Packed unit / Weight	0.407	Bare Unit / Weight (kg)	0.407
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

03/04/2026 00:00:07 – For information purposes only, can be subject to change