

## S/S BOX 20 X 10 X 5 CM

Ref. AI20105



02/04/2026 23:12:05 – For information purposes only, can be subject to change

## S/S BOX 20 X 10 X 5 CM

Ref. [AI20105](#)



**Convenient for storing surgical instruments**

### Description

- Sold with lid to protect instruments
- Decontamination with detergent and disinfectant -
- Sterilizable by autoclave (134°-18 minutes)
- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- Available in different sizes

### Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019505	HS Code	90189084
Warranty (months)	12	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	11	Bare Unit / Width (cm)	10.6
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	5.1
Packed unit / Length (cm)	21	Bare Unit / Length (cm)	20.5
Packed unit / Weight	0.397	Bare Unit / Weight (kg)	0.397
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 02/04/2026 23:12:05 – For information purposes only, can be subject to change