

COMPLETE ENT. SURGERY SET

Ref. IBOI301



Description

IA02402 Tongue depressor, stainless steel, curved, 135° angle, adult, 16 cm long IC08516 Metzemaum scissors, curved, 16 cm
IC09411 Iridectomy scissors, straight, 11.5 cm
IC11414 Mousses scissors, straight, 14 cm
IC39118 Gillies hook, 18 cm
IC46127 Yankauer/Ombredanne cannula, 27 cm
ID03600 Doyen mouth opener, paediatric
ID03602 Doyen mouth opener, adult
ID05403 Scalpel handle no. 3
IE02612 Farabeuf retractor, 12 cm x 12 mm, pair
IE32217 Adson elevator, 17 cm
AI15080 Beak cup, stainless steel, dia. 80 mm, 165 ml
AI35158 Stainless steel box, 35 x 18 x 8 cm
IP02014 Kocher forceps, straight, A/G, 14 cm
IP08113 Halstead forceps, 13 cm, curved, A/G
IP08313 Halstead forceps, 13 cm, curved, S/G
IP10612 Backhaus field forceps, 12 cm
IP11613 Chaput forceps, 13 cm
IP14514 Hartmann sliding forceps, for foreign bodies, 14 cm, jaws 12 x 3 mm
P16514 Kelly forceps, S/G, curved, 14 cm
IP21414 Murphy Pean forceps, 14 cm
IP24516 Tilley-Lubet Barbon Nasal Forceps, 16 cm, serrated jaws
IP26014 Dissection Forceps, A/G, 14 cm
IP26214 Dissection Forceps, S/G, 14 cm
IP49020 Mixer O'Shaughnessy Forceps, 90°, 20 cm
IP50115 *Crile-Wood Needle Holder, 15 cm, tungsten
IR07506 Farabeuf Rasp, curved, 15 cm x 6 mm
IS17110 Ear Syringe, 110 ml
IS23100 Vacher Nasal Speculum, with window, adult, 9 mm
IS23200 Vacher Nasal Speculum, with window, child, 6 mm

COMPLETE ENT. SURGERY SET

Ref. IBOI301



Features

CE Marking	CE	Medical device	Class I
Bar code	3661040082400	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	3.31
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

03/04/2026 00:04:22 – For information purposes only, can be subject to change