

EYE EXAMINATION SET

Ref. IBOI401



02/04/2026 23:25:59 – For information purposes only, can be subject to change

EYE EXAMINATION SET

Ref. IBOI401



Description

This box contains:
 Stainless steel box, 20 x 10 x 5 cm -IS
 Pley blepharostat, 7 cm
 Desmarres eyelid retractor, 10 mm
 Magnet + Handle
 Graefe dissection forceps, fine, with iris, 10 cm, straight, S/G
 Iridectomy scissors, straight, 10 cm -IS, Holtex
 Mousse scissors, straight, 14 cm -IS, Holtex
 Graefe hook for strabismus, foam
 Hunter tweezers, square jaws, 8 cm -IS
 Dissection forceps, S/G, fine, pointed, 14 cm -IS, Holtex

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040082646	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.77

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 02/04/2026 23:25:59 – For information purposes only, can be subject to change