

HOOK TO REMOVE IUD 25CM (DOWNS P-347)

Ref. [IC38201](#)



Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040011448	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.026
Bare Unit / Weight (kg)	0.024		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

03/04/2026 01:35:26 – For information purposes only, can be subject to change