

BECKMANN RETRACTOR (4X4 TEETH) 21 CM

Ref. IE13421



CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2 ANS
2
2 YEARS
GUARANTEE

STAINLESS
INOX
STEEL

03/04/2026 01:33:08 – For information purposes only, can be subject to change

BECKMANN RETRACTOR (4X4 TEETH) 21 CM



Ref. IE13421

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040016634	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	10
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	2.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	20.4
Packed unit / Weight	0.123	Bare Unit / Weight (kg)	0.123
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

03/04/2026 01:33:08 – For information purposes only, can be subject to change