


**BERGMANN Mallet (2 HEADS - 510 G) 23
CM S/S**

Ref. IM02500



CE

DISPOSITIF MEDICAL
CLASSE
Ir
MEDICAL DEVICE

GARANTIE
2
2 YEARS
GUARANTEE

STAINLESS
INOX
STEEL

03/04/2026 01:53:27 – For information purposes only, can be subject to change

BERGMANN MALLETT (2 HEADS - 510 G) 23 CM S/S



Ref. IM02500

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040020419	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	6	Bare Unit / Width (cm)	6
Packed unit / Height (cm)	3	Bare Unit / Height (cm)	3.7
Packed unit / Length (cm)	23	Bare Unit / Length (cm)	23.3
Packed unit / Weight	0.607	Bare Unit / Weight (kg)	0.607

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

03/04/2026 01:53:27 – For information purposes only, can be subject to change