

CHAPUT FORCEPS 16 CM

Ref. IP11616




CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2 ANS
2
2 YEARS
GUARANTEE

STAINLESS
INOX
STEEL



03/04/2026 00:23:08 – For information purposes only, can be subject to change

CHAPUT FORCEPS 16 CM

Ref. IP11616



Features

CE Marking	CE	Medical device	Class I
Bar code	3661040161068	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.028
Bare Unit / Weight (kg)	0.028		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--