

MINI BAG OUessant, DE BOISSY, WITH BLACK LEATHER POCKET

Ref. MDOUES02

De Boissy
— DEPUIS 1982



📅 03/04/2026 02:01:15 – For information purposes only, can be subject to change

MINI BAG OUessant, DE BOISSY, WITH BLACK LEATHER POCKET

Ref. MDOUES02



De Boissy
— DEPUIS 1982

The mini version of the OUessant Medical Bag

Description

Eco-responsible, customisable and 'super-organised' for a trendy look with all the practical, modern options.

All the fabrics, both exterior and interior, are Global Recycled Standard (GRS) certified, made from recycled materials and remain washable, hard-wearing and water-repellent.

Its deep blue, inspired by the shores of Ushant, is enhanced by black grained leather-effect finishes, giving it an elegant, professional look.

While De Boissy was keen to pay tribute to carers and take inspiration for its bags from the shape of the traditional doctors leather satchel, the look of the Ouessant bags contrasts sharply with the usual classicism. De Boissy has opted for elegance, trendiness and modernity. And the fact remains that they have been designed to be 'ultra practical', to be the best ally of healthcare professionals during their medical visits at home or in the community.

- The mini version of the Etrepat and Eiffel Medical Bags
- The main space inside also includes 6 pockets, a divider and 2 mobile compartments: 1 transparent and 1 isothermal, supplied with its reusable cold gel pouch.
- The belly pocket measures 31×15 cm
- Still on the outside: the 3 small zipped pockets, plus the one with elastic for the gel bottle, the detachable-customisable grained leather-look pouch, and the external USB port with integrated charging cable
- Same handles, shoulder strap, name tag and non-slip bottom as the large format
- Dimensions (W x H x D): 32 x 28 x 20 cm
- Weight: 1.6 kg
- Maximum recommended contents: 4 kg / 15 litres
- Sold empty

The De Boissy bag comes in a dust bag to protect it when stored. Grey with patterns, this accessory is made from recycle materials with a drawstring fastening.

MINI BAG OUESSANT, DE BOISSY, WITH BLACK LEATHER POCKET

Ref. MDOUES02



Features

| | | | |
|------------------------------------|---------------|---------------------------|----------|
| CE Marking | NONE | Medical device | No |
| Bar code | 3700446042424 | HS Code | 42022290 |
| Warranty (months) | 12 | Country | China |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 54 | Packed unit / Height (cm) | 42 |
| Packed unit / Length (cm) | 57 | Packed unit / Weight | 1.05 |
| Bare Unit / Weight (kg) | 1.05 | | |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |

03/04/2026 02:01:15 – For information purposes only, can be subject to change