

ADULT SILICON SPO2 SENSOR FOR VITALIMATE MINI

Ref. [MONI010A](#)

The silicon SpO2 sensor for Vitalimate MINI

Description

Silicon SpO2 sensor for Vitalimate MINI

- adult size
- length 1 m
- 6 pin sensor

Available in pediatric size

Spengler also offers all sizes rigid SpO2 sensors for the Vitalimate MINI as well as an ear clip sensor..

Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	6948309002544	HS Code	9018193090
Warranty (months)	6	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.055
Bare Unit / Weight (kg)	0.05		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	10	Maximum humidity	93

OTHER FEATURES

03/04/2026 00:24:45 – For information purposes only, can be subject to change