

**Ref. MONI010D**



**CE**  
**0123**

DISPOSITIF MEDICAL  
CLASSE  
**II<sub>b</sub>**  
MEDICAL DEVICE

## Description

03/04/2026 00:51:42 – For information purposes only, can be subject to change


 Ref. **MONI010D**

## Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	6948309002537	HS Code	9018193090
Warranty (months)	6	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.055
Bare Unit / Weight (kg)	0.05		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	55
Minimum humidity	10	Maximum humidity	93
<b>OTHER FEATURES</b>			

03/04/2026 00:51:42 – For information purposes only, can be subject to change