

S/S BOX 25 X 10 X 5 CM

Ref. AI25105



06/04/2026 21:39:55 – For information purposes only, can be subject to change

S/S BOX 25 X 10 X 5 CM

Ref. [AI25105](#)



Ideal for storing surgical instruments

Description

- Sold with lid to protect instruments
-
- Decontamination with detergent and disinfectant -
-
- Sterilizable by autoclave (134°-18 minutes)
-
- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
-
- Available in different sizes

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019543	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	10	Bare Unit / Width (cm)	10
Packed unit / Height (cm)	5	Bare Unit / Height (cm)	5
Packed unit / Length (cm)	25	Bare Unit / Length (cm)	25
Packed unit / Weight	0.563	Bare Unit / Weight (kg)	0.563

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 06/04/2026 21:39:55 – For information purposes only, can be subject to change