

S/S KIDNEY DISH 30 X 14 X 4,5 CM

Ref. AI28124



06/04/2026 22:11:36 – For information purposes only, can be subject to change

S/S KIDNEY DISH 30 X 14 X 4,5 CM

Ref. [AI28124](#)



Ideal for storing instruments or temporarily installing medical waste

Description

- Decontamination with detergent and disinfectant -
- Autoclavable (134°-18 minutes)- In the " ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- Reference 100% IN-STORE: this service instrumentation reference is guaranteed available immediately
- Available in different sizes

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019567	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	17	Bare Unit / Width (cm)	14
Packed unit / Height (cm)	7	Bare Unit / Height (cm)	4.5
Packed unit / Length (cm)	21	Bare Unit / Length (cm)	30
Packed unit / Weight	0.28	Bare Unit / Weight (kg)	0.246
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 06/04/2026 22:11:36 – For information purposes only, can be subject to change