

**S/S SERVING FORCEPS JAR WITH
COVER(DIAM)50X200 MM**

Ref. AI11723



06/04/2026 21:23:49 – For information purposes only, can be subject to change

S/S SERVING FORCEPS JAR WITH COVER(DIAM)50X200 MM



Ref. AI11723

Features

| | | | |
|------------------------------------|---------------|-------------------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3661040019321 | HS Code | 90189084 |
| Warranty (months) | 12 | Country | Pakistan |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 11 | Bare Unit / Width (cm) | 5 |
| Packed unit / Height (cm) | 11 | Bare Unit / Height (cm) | 20 |
| Packed unit / Length (cm) | 24 | Bare Unit / Length (cm) | 20 |
| Packed unit / Weight | 0.222 | Bare Unit / Weight (kg) | 0.222 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

06/04/2026 21:23:49 – For information purposes only, can be subject to change