

LID FOR BEDPAN (STAINLESS STEEL)
48X30CM

Ref. AI14831



06/04/2026 20:49:31 – For information purposes only, can be subject to change

LID FOR BEDPAN (STAINLESS STEEL) 48X30CM



Ref. AI14831

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040032917	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	54	Bare Unit / Width (cm)	30
Packed unit / Height (cm)	2	Packed unit / Length (cm)	39
Bare Unit / Length (cm)	48	Packed unit / Weight	0.416
Bare Unit / Weight (kg)	0.416		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

06/04/2026 20:49:31 – For information purposes only, can be subject to change