

S/S BOWL WITH SPOUT DIAM 80 MM - 150 ML

Ref. AI15080



06/04/2026 21:36:44 – For information purposes only, can be subject to change

S/S BOWL WITH SPOUT DIAM 80 MM - 150 ML



Ref. AI15080

Allows surgical instruments to be stored and transported

Description

Holtex instruments are all made of premium stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to ensure good durability and corrosion resistance.

The steels come into contact with human tissues, so they must comply with international standards and regulations. Holtex guarantees you medical devices that comply with Directive 93/42/EC and ISO 13485 (2016).

- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual

- Reference 100% IN-STORE: this service instrumentation reference is guaranteed to be available immediately.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019383	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	17	Bare Unit / Width (cm)	8
Packed unit / Height (cm)	7	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	21	Bare Unit / Length (cm)	8
Packed unit / Weight	0.056	Bare Unit / Weight (kg)	0.056

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 06/04/2026 21:36:44 – For information purposes only, can be subject to change