

**S/S BOX PERFORATED 20 X 10 X 3 CM
PERFORATED TOP / BOTTOM**

Ref. AI20103P



06/04/2026 21:25:52 – For information purposes only, can be subject to change

S/S BOX PERFORATED 20 X 10 X 3 CM PERFORATED TOP / BOTTOM



Ref. AI20103P

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040067735	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	13	Bare Unit / Width (cm)	13.6
Packed unit / Height (cm)	5	Bare Unit / Height (cm)	5.5
Packed unit / Length (cm)	28	Bare Unit / Length (cm)	28
Packed unit / Weight	0.325	Bare Unit / Weight (kg)	0.325
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

06/04/2026 21:25:52 – For information purposes only, can be subject to change