

**S/S BOX PERFORATED 20 X 10 X 6 CM  
PERFORATED TOP / BOTTOM**

Ref. AI20106P



06/04/2026 20:24:10 – For information purposes only, can be subject to change

# S/S BOX PERFORATED 20 X 10 X 6 CM PERFORATED TOP / BOTTOM



Ref. AI20106P

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040067759	HS Code	90189084
Warranty (months)	12	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	6	Bare Unit / Width (cm)	6.9
Packed unit / Height (cm)	3	Bare Unit / Height (cm)	3.2
Packed unit / Length (cm)	6	Bare Unit / Length (cm)	6.9
Packed unit / Weight	0.391	Bare Unit / Weight (kg)	0.391
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

06/04/2026 20:24:10 – For information purposes only, can be subject to change