

**STAINLESS STEEL BOX, 25 X 10 X 5 CM -
PERFORATED TOP/BOTTOM**

Ref. AI25105P



06/04/2026 19:30:20 – For information purposes only, can be subject to change

STAINLESS STEEL BOX, 25 X 10 X 5 CM - PERFORATED TOP/BOTTOM



Ref. [AI25105P](#)

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040067766	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	26	Bare Unit / Width (cm)	10
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	5
Packed unit / Length (cm)	12	Bare Unit / Length (cm)	25
Packed unit / Weight	6.7		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

06/04/2026 19:30:20 – For information purposes only, can be subject to change