

STAINLESS STEEL INSTRUMENT TRAY 40X30X4 CM

Ref. AI40304



06/04/2026 20:11:36 – For information purposes only, can be subject to change

STAINLESS STEEL INSTRUMENT TRAY 40X30X4 CM



Ref. AI40304

Ideal for storing small surgical instruments

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040067919	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	35	Bare Unit / Width (cm)	30
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	46	Bare Unit / Length (cm)	40
Packed unit / Weight	1	Bare Unit / Weight (kg)	0.9

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 06/04/2026 20:11:36 – For information purposes only, can be subject to change