

RETINAL DETACHMENT SURGERY SET

Ref. IBOI407



Description

Features

| | | | |
|------------------------------------|---------------|---------------------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3661040082707 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | Pakistan |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | 0.6 |
| Bare Unit / Weight (kg) | 0.6 | | |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

06/04/2026 21:29:21 – For information purposes only, can be subject to change