

# ROCHESTER-PEAN FORCEPS CURVED 18 CM

Ref. IP02318



CE

DISPOSITIF MEDICAL  
CLASSE  
**I**  
MEDICAL DEVICE

GARANTIE  
2 ANS  
**2**  
2 YEARS  
GUARANTEE

STAINLESS  
**INOX**  
STEEL

06/04/2026 21:15:02 – For information purposes only, can be subject to change

# ROCHESTER-PEAN FORCEPS CURVED 18 CM



Ref. IP02318

**Convenient to stop bleeding without damaging fragile tissues during an operation. Also used to perform care (dressings, etc.).**

## Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

100% IN-STORE reference: this reference of service instrumentation is guaranteed to be available immediately.

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040021652	HS Code	90189084
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	23	Bare Unit / Width (cm)	8.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.4
Packed unit / Length (cm)	9	Bare Unit / Length (cm)	18
Packed unit / Weight	0.06	Bare Unit / Weight (kg)	0.061
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 06/04/2026 21:15:02 – For information purposes only, can be subject to change