

CORWIN TONSIL FORCEPS

Ref. IC99300



CE



06/04/2026 19:53:09 – For information purposes only, can be subject to change

CORWIN TONSIL FORCEPS

Ref. IC99300



Features

| | | | |
|------------------------------------|---------------|-------------------------|----------|
| CE Marking | CE | Medical device | Class Ir |
| Bar code | 3661040013244 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | Pakistan |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Bare Unit / Width (cm) | 6.7 |
| Packed unit / Height (cm) | nc | Bare Unit / Height (cm) | 2.5 |
| Packed unit / Length (cm) | nc | Bare Unit / Length (cm) | 16.6 |
| Packed unit / Weight | 0.083 | Bare Unit / Weight (kg) | 0.083 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

06/04/2026 19:53:09 – For information purposes only, can be subject to change