

HARTMANN RETRACTOR STAINLESS STEEL WIRE 16 CM (PAIR)



Ref. IE03616

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040016375	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	5.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	16
Packed unit / Weight	0.1	Bare Unit / Weight (kg)	0.085
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 06/04/2026 22:14:28 – For information purposes only, can be subject to change