

BECKMANN RETRACTOR (4X4 DIENTES) 16 CM

Ref. IE13416



CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2 ANS
2
2 YEARS
GUARANTEE

STAINLESS
INOX
STEEL

06/04/2026 22:11:56 – For information purposes only, can be subject to change

BECKMANN RETRACTOR (4X4 DIENTES) 16 CM



Ref. IE13416

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040016627	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	8.3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	2
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	16.4
Packed unit / Weight	0.088	Bare Unit / Weight (kg)	0.088
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

06/04/2026 22:11:56 – For information purposes only, can be subject to change