

FINOCHIETTO LIGATURE CARRIER 24 CM

Ref. IP44524



Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040025582	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	8.6
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	24.7
Packed unit / Weight	0.088	Bare Unit / Weight (kg)	0.088
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

06/04/2026 19:47:48 – For information purposes only, can be subject to change