

# STAINLESS STEEL BED PAN SLIPPER TYPE 42X29 CM WITHOUT LID

Ref. AI14832



06/04/2026 19:50:59 – For information purposes only, can be subject to change

# STAINLESS STEEL BED PAN SLIPPER TYPE 42X29 CM WITHOUT LID



Ref. AI14832

## Description

- Maximum weight supported: 90 kg

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040115207	HS Code	90189084
Warranty (months)	12	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	30	Bare Unit / Width (cm)	29
Packed unit / Height (cm)	17	Packed unit / Length (cm)	56
Bare Unit / Length (cm)	42	Packed unit / Weight	1.08
Bare Unit / Weight (kg)	1.481		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

06/04/2026 19:50:59 – For information purposes only, can be subject to change