

S/S TRAY WITH COVER 280 X 180 X 30 MM

Ref. AI30214



06/04/2026 20:12:19 – For information purposes only, can be subject to change

S/S TRAY WITH COVER 280 X 180 X 30 MM

Ref. AI30214



Ideal for storing small surgical instruments

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040039022	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	21	Bare Unit / Width (cm)	21.5
Packed unit / Height (cm)	5	Bare Unit / Height (cm)	5
Packed unit / Length (cm)	30	Bare Unit / Length (cm)	30
Packed unit / Weight	0.95	Bare Unit / Weight (kg)	0.95
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 06/04/2026 20:12:19 – For information purposes only, can be subject to change