

**STAINLESS STEEL INSTRUMENT TRAY 46.5
X 31 X 5 CM**

Ref. AI46315



Ideal for storing small surgical instruments

06/04/2026 20:57:49 – For information purposes only, can be subject to change

STAINLESS STEEL INSTRUMENT TRAY 46.5 X 31 X 5 CM



Ref. AI46315

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040061412	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	33	Packed unit / Height (cm)	6
Packed unit / Length (cm)	48	Packed unit / Weight	1
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

06/04/2026 20:57:49 – For information purposes only, can be subject to change