


# CHAPUT FORCEPS 18 CM

Ref. IP11618



06/04/2026 19:24:37 – For information purposes only, can be subject to change

# CHAPUT FORCEPS 18 CM

Ref. IP11618



## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040161075	HS Code	90189084
Warranty (months)	24	Country	Pakistan

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.028
Bare Unit / Weight (kg)	0.028		

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--

06/04/2026 19:24:37 – For information purposes only, can be subject to change