

# LUER BONE RONGEUR MULTIACTION STRAIGHT 22 CM

Ref. IP48422



06/04/2026 20:03:54 – For information purposes only, can be subject to change

# LUER BONE RONGEUR MULTIACTION STRAIGHT 22 CM



Ref. IP48422

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040027098	HS Code	90189084
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	9.7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	22.7
Packed unit / Weight	0.364	Bare Unit / Weight (kg)	0.364
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

06/04/2026 20:03:54 – For information purposes only, can be subject to change