

PERCY AMPUTATION RETRACTOR (DIAM 17 CM) REMOVABLE HANDLES

Ref. IR05100



03/04/2026 02:05:13 – For information purposes only, can be subject to change

PERCY AMPUTATION RETRACTOR (DIAM 17 CM) REMOVABLE HANDLES



Ref. IR05100

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040029030	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	19	Bare Unit / Width (cm)	19
Packed unit / Height (cm)	2	Bare Unit / Height (cm)	2
Packed unit / Length (cm)	19	Bare Unit / Length (cm)	19
Packed unit / Weight	0.6	Bare Unit / Weight (kg)	0.599
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

03/04/2026 02:05:13 – For information purposes only, can be subject to change