

HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°30

Ref. IS01630



Exemple de bougies Hégar simple, inox.
30 tailles disponibles, nous consulter



06/04/2026 21:32:10 – For information purposes only, can be subject to change

HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°30



Ref. IS01630

Description

Hégar Single Ended Dilator Stainless Steel

30 Sizes available, contact us

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040030005	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	2.6
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	2.6
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	19
Packed unit / Weight	0.227	Bare Unit / Weight (kg)	0.227

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

06/04/2026 21:32:10 – For information purposes only, can be subject to change