

METAL FEMALE CATHETER 15 CM

Ref. IS03402



Description

4mm diameter

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040030128	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	1.1
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.7
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	16.2
Packed unit / Weight	0.012	Bare Unit / Weight (kg)	0.012

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

06/04/2026 20:02:15 – For information purposes only, can be subject to change