

**PINARD OBSTETRICAL STETHOSCOPE,  
ALUMINIUM**

Ref. IS07601



CE



06/04/2026 22:16:32 – For information purposes only, can be subject to change

# PINARD OBSTETRICAL STETHOSCOPE, ALUMINIUM

Ref. IS07601



**Stethoscope foetal Pinard used to listen to foetal heart sounds**

## Description

Small conical instrument, of which the top allows to the practitioner to place the ear, while the flat and circular basis is placed on the abdomen of the pregnant woman.  
Designed for frequent and disinfection with hospital-grade products.

- Single chestpiece
- Lightweight
- Unbreakable aluminum
- Easy to handle
- Easy to clean

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040030548	HS Code	90189084
Warranty (months)	24	Country	Pakistan

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	13	Bare Unit / Width (cm)	7
Packed unit / Height (cm)	7	Bare Unit / Height (cm)	7
Packed unit / Length (cm)	23	Bare Unit / Length (cm)	14
Packed unit / Weight	0.034	Bare Unit / Weight (kg)	0.034

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 06/04/2026 22:16:32 – For information purposes only, can be subject to change