

OCHSNER TROCART 14 FRENCH GAUGE (4.66 MM)

Ref. IT11614



Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040031729	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	7	Packed unit / Height (cm)	1
Packed unit / Length (cm)	21	Packed unit / Weight	0.028
Bare Unit / Weight (kg)	0.028		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

06/04/2026 20:49:06 – For information purposes only, can be subject to change