

# PRINTING PAPER 50MM\*20M (NO LOGO) FOR NC5 (BOX OF 10 ROLLS)

Ref. BCIMPR02

High-quality ECG paper rolls for precise and continuous printing, compatible with the Comen NC5.

## Description

High-quality ECG paper rolls for precise and continuous printing, compatible with the Comen NC5.

## Features

CE Marking	NONE	Medical device	No
Bar code	3661040219394	HS Code	482340000
Country	China		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 06/04/2026 22:18:03 – For information purposes only, can be subject to change