



**Ref. MONI011E**



CE



06/04/2026 19:35:13 – For information purposes only, can be subject to change

=====

**Ref. MONI011E**

=====

## Description

=====

## Features

CE Marking	CE	Medical device	Class I
Bar code	6948309002582	HS Code	9018193090
Warranty (months)	6	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	2
Packed unit / Length (cm)	18	Packed unit / Weight	0.15
Bare Unit / Weight (kg)	0.1		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	55
Minimum humidity	10	Maximum humidity	93
<b>OTHER FEATURES</b>			

 06/04/2026 19:35:13 – For information purposes only, can be subject to change