

S/S BOX 20 X 10 X 6 CM

Ref. AD20106



07/04/2026 02:17:45 – For information purposes only, can be subject to change

S/S BOX 20 X 10 X 6 CM

Ref. [AD20106](#)



Convenient for storing surgical instruments

Description

- Sold with lid to protect instruments
- Decontamination with detergent and disinfectant -
- Sterilizable by autoclave (134°-18 minutes)
- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- Available in different sizes

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040013657	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	11	Bare Unit / Width (cm)	10.7
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	6
Packed unit / Length (cm)	21	Bare Unit / Length (cm)	20.5
Packed unit / Weight	0.451	Bare Unit / Weight (kg)	0.451

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 07/04/2026 02:17:45 – For information purposes only, can be subject to change