

S/S BOX 28 X 14 X 6 CM

Ref. AI28146



07/04/2026 02:37:31 – For information purposes only, can be subject to change

S/S BOX 28 X 14 X 6 CM

Ref. [AI28146](#)



Convenient for storing surgical instruments

Description

- Sold with lid to protect instruments
- Decontamination with detergent and disinfectant -
- Sterilizable by autoclave (134°-18 minutes)- In the " ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- Reference 100% IN-STORE: this service instrumentation reference is guaranteed available immediately
- Available in different sizes

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019574	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	14	Bare Unit / Width (cm)	14
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	6
Packed unit / Length (cm)	28	Bare Unit / Length (cm)	28
Packed unit / Weight	0.872	Bare Unit / Weight (kg)	0.872

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 07/04/2026 02:37:31 – For information purposes only, can be subject to change