

S/S BOX 35 X 18 X 8 CM

Ref. AI35158



06/04/2026 23:41:37 – For information purposes only, can be subject to change

S/S BOX 35 X 18 X 8 CM

Ref. [AI35158](#)



Convenient for storing surgical instruments

Description

- Sold with lid to protect instruments
- Decontamination with detergent and disinfectant -
- Sterilizable by autoclave (134°-18 minutes)
- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- Available in different sizes

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019659	HS Code	90189084
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	18	Bare Unit / Width (cm)	18
Packed unit / Height (cm)	8	Bare Unit / Height (cm)	8
Packed unit / Length (cm)	35	Bare Unit / Length (cm)	35
Packed unit / Weight	1.25	Bare Unit / Weight (kg)	1.25
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

06/04/2026 23:41:37 – For information purposes only, can be subject to change