

## FHR1 PROBE/WATER PROOF

Ref. BCTOCO13

COMEN



CE  
0482



### C11 cardiocograph probe to monitor the baby's heart rate

#### Description

Very sensitive probes (9 crystals)  
Waterproof probes (IP67), ideal for childbirth

06/04/2026 23:49:46 – For information purposes only, can be subject to change

# FHR1 PROBE/WATER PROOF

Ref. BCTOCO13

# COMEN

## Features

CE Marking	CE0482	Medical device	Class IIa
Bar code	3661040219240	HS Code	9018193090
Warranty (months)	6	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

06/04/2026 23:49:46 – For information purposes only, can be subject to change