

FHR2 PROBE/ WATER PROOF

Ref. BCTOCO14

COMEN



CE
0482



C11 cardiocograph probe to monitor the baby's crdiac rhythm.

Description

Very sensitive probes (9 crystals)
Waterproof probes (IP67), ideal for childbirth.

📅 07/04/2026 02:08:59 – For information purposes only, can be subject to change

FHR2 PROBE/ WATER PROOF

COMEN

Ref. BCTOCO14

Features

CE Marking	CE0482	Medical device	Class IIa
Bar code	3661040219257	HS Code	9018193090
Warranty (months)	6	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

07/04/2026 02:08:59 – For information purposes only, can be subject to change