

TOCO PROBE / WATER PROOF

Ref. BCTOCO15

COMEN



CE
0482



TOCO probe of the cardiocotograph to monitor the contractions of the mother-to-be

Description

IP67 waterproof

 07/04/2026 00:05:08 – For information purposes only, can be subject to change

TOCO PROBE / WATER PROOF

Ref. BCTOCO15

COMEN

Features

CE Marking	CE0482	Medical device	Class IIa
Bar code	3661040219264	HS Code	9018193010
Warranty (months)	6	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

07/04/2026 00:05:08 – For information purposes only, can be subject to change