

LARYNGAL MIRROR N°9 28 MM

Ref. IM03109



07/04/2026 01:00:44 – For information purposes only, can be subject to change

LARYNGAL MIRROR N°9 28 MM

Ref. IM03109



Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

Features

| | | | |
|-------------------|---------------|----------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3661040020501 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | Pakistan |

DIMENSIONS AND WEIGHT

| | | | |
|---------------------------|------|-------------------------|------|
| Packed unit / Width (cm) | nc | Bare Unit / Width (cm) | 2.8 |
| Packed unit / Height (cm) | nc | Bare Unit / Height (cm) | 1.3 |
| Packed unit / Length (cm) | nc | Bare Unit / Length (cm) | 14.6 |
| Packed unit / Weight | 0.01 | Bare Unit / Weight (kg) | 0.01 |

SPECIFIC STORAGE CONDITIONS

| | | | |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |

OTHER FEATURES

| | |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|

07/04/2026 01:00:44 – For information purposes only, can be subject to change